

PRODUCT CATALOG

Emsculpt And Emsella Series	01 - 06
Cryolipolysis Series	07 - 26
Spa Series	27 - 48
Diode Laser Hair Removal Series	49 - 60
Pico Laser Series	61 - 66
IPL & ND YAG Laser Tattoo Removal Series	67 - 70
Wrinkle Removal Series	71 - 76
Velashape Body Slimming Machine	77 - 82
Lipolaser Series	83 - 96
Pressotherapy Series	97 - 106
Shock Wave Therapy Series	107 - 108
Multi-function Weight Loss Series	109 - 116
Laser Hair Regrowth Series	117 - 120
Icevind Skin Cooling machine	121 - 122
Mini Home Use Beauty Machine	123 - 124
Derma Roller Series & Tattoo Gun	125 - 130

EMSCULPT AND EMSELLA SERIES

WORKING PRINCIPLE


Using high-intensity focused electromagnetic field technology → through stimulation of nerves

- → continuous expansion and contraction of autologous muscles → extreme training
- → deep remodeling of muscles→ growth of myofibrils (muscle enlargement)
- → production of new collagen chains and muscle fibers (muscles) Hyperplasia)
- → to train and increase muscle density and volume



Using focused magnetic resonance (HIEMT) technology→100% limit muscle contraction


- →can trigger a large amount of lipolysis →fatty acid is broken down from triglyceric acid
- →accumulate in a large amount in fat cells → make fatty acid concentration too high
- →adipocyte apoptosis→ Excreted from the body by normal meta-bolism within 2-4 weeks





High-Intensity Focused Electromagnetic handle 16*23CM

O SCOPE OF APPLICATION







Hip

Thigh

Waist

FUNCTIONS

- 1). Exercise abdominal muscles Shape the vest line
- 2. Exercise buttocks muscles Shape the peach hip
- 3. Exercise the oblique muscles Shape the Mermaid Line
- Solve postpartum Rectus abdominis separation
- (5). Treatment of waist muscle strain

(1) COMPARISON OF TREATMENT EFFECTS









BEFORE



AFTER

BEFORE

<01>

High Frequency Body Muscle Stimulation Electromagnetic Vibration HIEMTSURE Body Sculpting Machine

Shape the perfect figure

PARAMETER

Model: MS-68

Voltage: AC110V-220V

Output intensity: 7 Tesla

Output Power: ≤2800W

Output Frequency: Max 100Hz

Fuse: 10AV15A

Machine Size: 111cm*65cm*42cm

Packing Size: 116cm*55cm*62cm

Gross Weight: 45KG



PARAMETER

Model	EMS-27
Voltage	AC110V-220V
Output intensity	7 Tesla
Output power	≤2800W
Output frequency	1-100Hz
Elastic bands	3
Fuse	10A/15A
Machine size	48.2cm*42.3cm*106cm
Packing size	58cm*73cm*117cm



PARAMETER

Model	MS-49
Туре	portable
Function	muscle building,hip lift,body sculpting
Treatment parts	abdomen, buttock, triceps
Voltage	AC110V-220V
Output intensity	7 Tesla
Output Power	≤2000W
Output Frequency	1-50Hz
Fuse	10A/15A
Machine Size	55cm*45cm*52cm
Packing Size	70cm*59cm*68cm
Gross Weight	40KG





WORKING PRINCIPLE



Pelvic floor muscles insufficiently support pelvic organs and affect bladder control.

EMSVela effectively stimulates pelvic floor muscles.

Stimulation leads to regained control over pelvic floor muscles and bladder.

S FUNCTIONS

- 1). Improve Frequent urination 2. Improve Incontinence
- 3. Improve prostate 4. Enhanced sexual function
- Promote postpartum recovery
- 6. Restore vaginal muscle tension and improve vaginal health
- Strengthen the pubococcygeal muscle

PARAMETER

Model	MS-35S
Output Intensity	2.5T
Output Power	300w-4500w
Output Frequency	1-60HZ
Fuse	15A
Treatment time	recommend 30 minutes
Host size	111*38*58cm
Chair Size	48*48*10cm
Host weight	45kg
Chair weight	3kg

